

Tailor-made primary packaging for pedicle screws

Developed in cooperation with medical technology companies and clinical specialists

- high quality and stable single use packaging for polyaxial and monoaxial pedicle screws
- optimum security and protection for product and user
- packaging tube for safe removal of the screw using the attaching tool
- can be used as primary package
- fixing of the pedicle screw by customer-/ product-specific clamping inside the tube
- speedy and efficient packing of products
- safe transport

Raw Material

- PETG

Packing

- tube and cap are packed separately (2 poly bags individually heat sealed)

Production environment

- ISO class 8 clean room

rose plastic medical packaging GmbH

Rupolzer Str. 30
D-88138 Hergensweiler/Lindau
Tel. ++49/8388/92345-0
Fax ++49/8388/92345-150
Info@rose-medipack.de
www.rose-medipack.de



customer- / product-specific clamping inside the tube

Article	Article description	Inside dimensions Ø (mm)	Usable length (mm)
Tube 1	short tube with product-specific clamping	up to approx. 17 mm	up to approx. 60 mm
Tube 2	long tube with product-specific clamping	up to approx. 17 mm	up to approx. 110 mm
Cap 1	cap without holding device for set screw		